EVIDENCIA vs. EXPERIENCIA

**TEG<sup>®</sup> 6s**

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### Índice

Introducción	3
<b>ARTÍCULOS SELECCIONADOS</b>	
Efficacy of point-of-care thromboelastography 6s to evaluate platelet function in a patient with pseudothrombocytopenia undergoing cardiopulmonary bypass: a case report	5
TEG6s Platelet Mapping assay for the estimation of plasma fibrinogen concentration during cardiovascular surgery: a single-center prospective observational study	7
Detection of hypofibrinogenemia during cardiac surgery: a comparison of resonance-based thrombelastography with the traditional Clauss method	9
Assessing TEG6s reliability between devices and across multiple time points: A prospective thromboelastography validation study	11
Coagulation ability when separating from cardiopulmonary bypass with and without fresh frozen plasma: a pilot study	13
Bibliografía	15

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### National Institute for Health and Care Excellence (NICE)

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## Efficacy of point-of-care thromboelastography 6s to evaluate platelet function in a patient with pseudothrombocytopenia undergoing cardiopulmonary bypass: a case report

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
Platelet count ( $\times 10^9/L$ )	Preoperative	During cardiopulmonary bypass	After protamine administration
EDTA	4.9	3.5	3.7
ACD	8.8	–	–
Rapid test from blood sampling (EDTA)	11.7	–	—
Rapid test from blood sampling (heparin)	10.4	–	
Fibrinogen (mg/dl)	369	211	226
Hemoglobin (g/dl)	10.9	6.7–8.2	8.4
Activated clotting time	103	427–700	111
Maximum amplitude of rapid TEG (mm, normal range: 52-70)	66.5	-	57.6
Amplitude at 10 min of rapid TEG (mm, normal range: 44-67)	63.1		48.8
Alpha angle of rapid TEG (degrees, normal range: 60-78)	76.2	-	69.7
Maximum amplitude of citrated functional fibrinogen (mm, normal range: 15-32)	27.6	-	20.0

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## Referencia

Chiba Y, Otsuka Y, Lefor AK, Sanui M. Efficacy of point-of-care thromboelastography 6s to evaluate platelet function in a patient with pseudothrombocytopenia undergoing cardiopulmonary bypass: a case report. *JA Clin Rep.* 2022;8(1):8.

## **TEG6s** Platelet Mapping assay for the estimation of plasma fibrinogen concentration during cardiovascular surgery: a single-center prospective observational study

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
## Referencia

Yamamoto Y, Sato Y, Takahashi M, Yamamoto H, Echizen M, Uchida T. TEG6s Platelet Mapping assay for the estimation of plasma fibrinogen concentration during cardiovascular surgery: a single-center prospective observational study. *J Anesth.* 2022;36(1):79-88.

# Detection of hypofibrinogenemia during cardiac surgery: a comparison of resonance-based thrombelastography with the traditional Clauss method

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

## Referencia

Magunia H, Azizy E, Krautter L, Rosenberger P, Straub A. Detection of hypofibrinogenemia during cardiac surgery: a comparison of resonance-based thrombelastography with the traditional Clauss method. Blood Coag Fibrinolysis. 2020;31(8):551-557.

# Assessing TEG6s reliability between devices and across multiple time points: A prospective thromboelastography validation study

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
## Referencia

Lloyd-Donald P, Churilov L, Cheong B, Bellomo R, McCall PR, Mårtensson J, Glassford N, Weinberg L. Assessing TEG6S reliability between devices and across multiple time points: A prospective thromboelastography validation study. *Sci Rep.* 2020;10(1):7045.

# Coagulation ability when separating from cardiopulmonary bypass with and without fresh frozen plasma: a pilot study.

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
## Referencia

Tamura T, Yokota S, Ito T, Ando M, Kubo Y, Waters JH, Nishiwaki K. Coagulation ability when separating from cardiopulmonary bypass with and without fresh frozen plasma: a pilot study. *Gen Thorac Cardiovasc Surg.* 2020;68(12):1361-1368.

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1. Varghese S, Unni M, Mukundan N, Ramji RL. Platelet functions in cardiopulmonary bypass surgery. *Med J Armed Forces India* 2005;61:316-321.
2. Dias JD, Lopez-Espina CG, Bleden K, Gurbel P, Hartmann J, Achneck HE. TEG®6s system measures the contributions of both platelet count and platelet function to clot formation at the site-of-care. *Platelets*. 2020;31(7):932-938.
3. Lloyd-Donald P, Vasudevan A, Angus P, Gow P, Mårtensson J, Glassford N, et al. Coagulation in acutely ill patients with severe chronic liver disease: Insights from thromboelastography. *J Crit Care*. 2017;38:215-224.
4. Fleming K, Redfern RE, March RL, Bobulski N, Kuehne M, Chen JT, et al. TEG-Directed Transfusion in Complex Cardiac Surgery: Impact on Blood Product Usage. *J Extra Corpor Technol*. 2017;49(4):283-290.
5. Nakayama Y, Nakajima Y, Tanaka KA, Sessler DI, Maeda S, Iida J, et al. Thromboelastometry-guided intraoperative haemostatic management reduces bleeding and red cell transfusion after paediatric cardiac surgery. *Br J Anaesth*. 2015;114:91-102.
6. Weber CF, Gorlinger K, Meiningner D, Herrmann E, Bingold T, Moritz A, et al. Point-of-care testing: a prospective, randomized clinical trial of efficacy in coagulopathic cardiac surgery patients. *Anesthesiology*. 2012;117:531-47.
7. Wikkelsø A, Wetterslev J, Møller AM, Afshari A. Thromboelastography (TEG) or thromboelastometry (ROTEM) to monitor haemostatic treatment versus usual care in adults or children with bleeding. *Cochrane Database Syst Rev*. 2016;(8):CD007871.
8. Korpálová B, Samoš M, Bolek T, Škorňová I, Kovář F, Kubisz P, et al. Role of Thromboelastography and Rotational Thromboelastometry in the Management of Cardiovascular Diseases. *Clin Appl Thromb Hemost*. 2018;24(8):1199-1207.
9. National Institute for Health and Care Excellence (NICE) Diagnostics Guidance 13: Detecting, managing and monitoring haemostasis: Viscoelastometric point-of-care testing (ROTEM, TEG and Sonoclot systems), 2014. Disponible en <https://www.nice.org.uk/guidance/dg13>. Fecha de acceso: noviembre 2022.
10. Pagano D, Milojevic M, Meesters MI, Benedetto U, Bolliger D, von Heymann C, et al. EACTS/EACTA Guidelines on patient blood management for adult cardiac surgery. *Eur J Cardiothorac Surg* 2018;53:79-111.
11. Raphael J, Mazer CD, Subramani S, Schroeder A, Abdalla M, Ferreira R, et al. Society of Cardiovascular Anesthesiologists Clinical Practice Improvement Advisory for Management of Perioperative Bleeding and Hemostasis in Cardiac Surgery Patients. *Anesth Analg*. 2019;129(5):1209-1221.
12. Bleden KP, Chaudhary R, Mohammed N, Muresan AA, Lopez-Espina CG, et al. Determination of non-Vitamin K oral anticoagulant (NOAC) effects using a new-generation thromboelastography TEG 6s system. *J Thromb Thrombolysis*. 2017;43(4):437-445.
13. Dias JD, Haney EI, Mathew BA, Lopez-Espina CG, Orr AW, Popovsky MA. New-Generation Thromboelastography: Comprehensive Evaluation of Citrated and Heparinized Blood Sample Storage Effect on Clot-Forming Variables. *Arch Pathol Lab Med*. 2017;141:569-577.
14. Ziegler B, Voelckel W, Zipperle J, Grottko O, Schöchl H. Comparison between the new fully automated viscoelastic coagulation analysers TEG 6s and ROTEM Sigma in trauma patients: A prospective observational study. *Eur J Anaesthesiol*. 2019;36(11):834-842.
15. Wong Q, Byrne KP, Robinson SC. Clinical agreement and interchangeability of TEG5000 and TEG6s during cardiac surgery. *Anaesth Intensive Care*. 2020;48(1):43-52.
16. Erdoes G, Schloer H, Eberle B, Nagler M. Next generation viscoelasticity assays in cardiothoracic surgery: Feasibility of the TEG6s system. *PLoS One*. 2018;13(12):e0209360.

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© 2015 Haemonetics Corporation. Haemonetics, Haemonetics The Blood Management Company, TEG y RapidTEG son marcas comerciales registradas comerciales de Haemonetics Corporation de los EE.UU. y/o en otros países.
Haemonetics Italia S.r.l.
Todos los derechos reservados. 11.2015 EE.UU. COL-COPIA-000923-ES(AB)