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1. Wiviott SD, Braunwald E, McCabe CH et al. Prasugrel versus clopidogrel in patients with acute coronary syndromes. *N Engl J Med* 2007; 357:2001–2015.
2. Wallentin L, Becker RC, Budaj A et al. Ticagrelor versus clopidogrel in patients with acute coronary syndromes. *N Engl J Med* 2009; 361: 1045–1057.
3. Curzen N. Antiplatelet therapy in acute coronary syndromes: beyond aspirin and clopidogrel. *Heart* 2012; 98: 1617-1619.
4. Sibbing D, Aradi D, Jacobshagen C et al. Guided de-escalation of antiplatelet treatment in patients with acute coronary syndrome undergoing percutaneous coronary intervention (TROPICAL-ACS): a randomised, open-label, multicentre trial. *Lancet* 2017; 390(10104):1747-1757.
5. Sousa-Uva M, Neumann FJ, Ahlsson A et al. 2018 ESC/EACTS Guidelines on myocardial revascularization. The Task Force on myocardial revascularization of the European Society of Cardiology (ESC) and European Association for Cardio-Thoracic Surgery (EACTS). *Eur J Cardiothorac Surg* 2019; 55:4-90. [Advance Access Publication]
6. Levine GN, Blankenship JC, Bailey SR et al. 2011 ACCF/AHA/SCAI Guideline for Percutaneous Coronary Intervention. *J Am Coll Cardiol* 2011; 58(24): e44–122.
7. Windecker S, Kohl P, Alfonso F et al. 2014 ESC/EACTS Guidelines on myocardial revascularization. The Task Force on myocardial revascularization of the European Society of Cardiology (ESC) and European Association for Cardio-Thoracic Surgery (EACTS). *Eur Heart J* 2014; 35:2541-2619.
8. Tantry US, Bonello L, Aradi D et al. Consensus and update on the definition of on-treatment platelet reactivity to adenosine diphosphate associated with ischemia and bleeding. *J Am Coll Cardiol* 2013; 62(24):2261-73.
9. Aradi D, Kirtane A, Bonello L et al. Bleeding and stent thrombosis on P2Y12-inhibitors: collaborative analysis on the role of platelet reactivity for risk stratification after percutaneous coronary intervention. *Eur Heart J* 2015; 36(27):1762-71.
10. Agarwal S, Johnson RI, Shaw M. Preoperative point-of-care platelet function testing in cardiac surgery. *J Cardiothorac Vasc Anesth* 2015; 29(2):333-341.
11. Kasivisvanathan R, Abbassi-Ghadi N, Kumar S et al. Risk of bleeding and adverse outcomes predicted by thromboelastography PlateletMapping in patients taking clopidogrel within 7 days of non-cardiac surgery. *Br J Surg* 2014; 101(11):1383-90.
12. Mahla E, Suarez TA, Bliden KP et al. Platelet function measurement-based strategy to reduce bleeding and waiting time in clopidogrel-treated patients undergoing coronary artery bypass graft surgery: the timing based on platelet function strategy to reduce clopidogrel-associated bleeding related to CABG (TARGET-CABG) study. *Circ Cardiovasc Interv* 2012; 5(2):261-9.
13. Sivapalan P, Back AC, Ostrowski SR et al. Transfusion requirements in elective cardiopulmonary bypass surgery patients: predictive value of Multiplate and Thromboelastography (TEG) PlateletMapping Assay. *Scand J Clin Lab Invest* 2017; 77(5):345-51.
14. Qi X, Zhao Y, Li K et al. Evaluating and monitoring the efficacy of recombinant activated factor VIIa in patients with haemophilia and inhibitors. *Blood Coagul Fibrinolysis* 2014; 25:754-60.
15. National Institute for Health and Care Excellence. Detecting, managing and monitoring haemostasis: viscoelastometric point-of-care testing (ROTEM, TEG and Sonoclot systems). NICE guidance 2014.
16. ASA Guidelines 2015: Practice guidelines for perioperative blood management: An updated report by the American Society of Anesthesiologists Task Force on perioperative blood management. *Anesthesiology* 2015; 122(2):241–275.
17. Rouillet S, de Maistre E, Ickx B et al. Position of the French Working Group on Perioperative Haemostasis (GIHP) on viscoelastic tests: What role for which indication in bleeding situations? *Anaesth Crit Care Med* 2018; S2352-5568(17)30283-7.
18. Patient blood management guidelines. Effect of perioperative strategies that minimize blood loss: Point-of-care testing. National Blood Authorities Australia.
19. Pagano D, Milojevic M, Meesters MI et al. 2017 EACTS/EACTA Guidelines on patient blood management for adult cardiac surgery. *Eur J Cardio Thorac Surg* 2018; 53:79–111.
20. Kozek-Langenecker SA, Afshari A, Albaladejo P et al. Management of severe perioperative bleeding: Guidelines from the European Society of Anaesthesiology. *Eur J Anaesthesiol* 2017; 34(6):332-395.
21. OAA/AAGBI Guidelines for Obstetric Anaesthetic Services 2013.
22. Collins P, Abdul-Kadir R, Thachil J, et al. Management of coagulopathy associated with postpartum hemorrhage: guidance from the SSC of the ISTH. *J Thromb Haemost* 2016; 14(1):205–210.
23. Rossaint R, Bouillon B, Cerny V et al. The European guideline on management of major bleeding and coagulopathy following trauma: fourth edition. *Critical Care* 2016; 20:100.

24. ACS TQIP Best Practice Guidelines 2012.
25. Gurbel PA, Bliden KP, Navickas IA et al. Adenosine diphosphate–induced platelet-fibrin clot strength: A new thrombelastographic indicator of long-term poststenting ischemic events. *Am Heart J* 2010; 160(2):346-54.
26. Sambu N, Radhakrishnan A, Dent H et al. Personalised antiplatelet therapy in stent thrombosis: observations from the Clopidogrel Resistance in Stent Thrombosis (CREST) registry. *Heart* 2012; 98(9):706-11.
27. Tang YD, Wang W, Yang M et al. Randomized comparisons of double-dose Clopidogrel or adjunctive Cilostazol versus standard Dual Antiplatelet in patients with high posttreatment platelet reactivity. Results of the CREATIVE trial. *Circulation* 2018; 137(21):2231-2245.
28. Mahla E, Tantry US, Prüller F et al. Is there a role for preoperative Platelet Function Testing in patients undergoing cardiac surgery during Antiplatelet Therapy? *Circulation* 2018; 138:2145–2159.
29. Kreuzt RP, Schmeisser G, Schaffter A et al. Prediction of ischemic events after Percutaneous Coronary Intervention: Thromboelastography profiles and Factor XIIIa activity. *TH Open* 2018; 2:e173-e181.
30. Cotton BA, Faz G, Hatch QM et al. Rapid thrombelastography delivers real-time results that predict transfusion within 1 hour of admission. *J Trauma* 2011; 71:407-14.
31. Introducing TEG®6s, 2015. Haemonetics Corp whitepaper.
32. Gurbel PA, Bliden KP, Tantry US et al. First report of the point-of-care TEG: A technical validation study of the TEG-6S system. *Platelets* 2016; 27(7):642-649.
33. Olechowski B, Dalton RT, Khanna V et al. Detection of individual responses to clopidogrel: Validation of a novel, rapid analysis using thrombelastography 6s. *Cardiovasc Ther* 2018; 36(4):e12433.
34. Dias JD, Hartmann J, Achneck HE. In vitro comparison of three whole blood platelet function analyzers: VerifyNow, Multiplate and PlateletMapping. ISTH 2019, Melbourne on July 6-10, 2019.

**Haemonetics S.A.**  
**Signy Centre, Rue des Fléchères 6**  
**1274 Signy-Centre**  
**Switzerland**

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© 2019 Haemonetics Corporation. Todos los derechos reservados. El contenido de este artículo técnico es propiedad de Haemonetics Corporation o sus filiales. Haemonetics y TEG son marcas comerciales o marcas registradas de Haemonetics Corporation en EE. UU., otros países o ambos. PlateletMapping es una marca de Cora Healthcare, Inc. 04.2019 Suiza. COL-PP-000252-IE(AA)